

The marker gene determines the preprocessing branch

Do not route fungal ITS and eukaryotic 18S through one bacterial 16S template.

Fungal ITS1 / ITS2

Locate variable spacer

ITSxpress · ITS1/ITS2 · taxa



Denoise variable lengths

DADA2 · trunc-len 0 if quality permits



Classify fungal ASVs

UNITE v10.0 · matching scope/date



Audit without forced tree

Length + chimera + non-target checks

Eukaryotic 18S SSU

Remove documented primers

Cutadapt · no ITSxpress



Denoise one fixed region

DADA2 · quality/overlap evidence



Classify eukaryotic ASVs

PR2 v5.1.1 · nine ranks



Filter after inspection

Plastid + bacteria + off-target audit